


# ProLink SV™

## Balloon Expandable Coronary Small Vessel Stent System



0.070 mm  
Articulation Width

0.090 mm  
Strut Thickness

0.30 - 0.35 bar  
Radial Strength

0.090 mm  
Strut Width

Manufactured By:  
 Vascular Concepts Limited  
 2 Deans Hall Business Park  
 Oak Road, Little Maplestead  
 Halstead, Essex C09 2RT  
 United Kingdom  
 Tel: 01787 477199  
 Fax: 01787 477510



Vascular Concepts Limited  
 No. 19, SV Complex (1st Floor)  
 Bellary Road, Hebbal  
 Bangalore 560 024  
 INDIA  
 Tel: (080) 23438146-49  
 Fax: (080) 23439205

## *ProLink SV Coronary Stent System Technical Data*

BALLOON [mm]	2.00	2.25	2.50
PRESSURE [Mpa/bar]			
0.6/6	2.00	2.25	2.50
0.7/7	2.01	2.26	2.51
0.8/8	2.02	2.27	2.52
0.9/9	2.04	2.29	2.54
1.0/10	2.05	2.31	2.56
1.1/11	2.07	2.33	2.59
1.2/12	2.09	2.35	2.61
1.3/13	2.11	2.37	2.64
1.4/14	2.13	2.39	2.66
1.5/15	2.14	2.41	2.69
1.6/16	2.16	2.43	2.72
1.7/17	2.18	2.45	2.75
1.8/18	2.20	2.48	2.79
1.9/19	2.22	2.50	2.81
2.0/20	2.24	2.53	2.83

NOMINAL PRESSURE AT 6 BAR (0.6MPa)

**RATED BURST PRESSURE**

## *ProLink SV Coronary Stent System Ordering Information*

LENGTH	BALLOON DIAMETER		
	2.0 mm	2.25 mm	2.50 mm
10 mm	PLSV2010	PLSV22510	PLSV2510
15 mm	PLSV2015	PLSV22515	PLSV2515
20 mm	PLSV2020	PLSV22520	PLSV2520
25 mm	PLSV2025	PLSV22525	PLSV2525
30 mm	PLSV2030	PLSV22530	PLSV2530